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	Výrobca: Famidoc Technology Co., Ltd. Room 401, 501, 601, 901, Building 6, Zhongtang Tian'an Cyber Park, No. 89, Jinyuan Road, Zhongtang Town, Dongguan City, Guangdong Province, 523220, P.R.China
---	--

	Zástupca pre EÚ: Názov: Shanghai International Holding Corp. GmbH (Evropa) Adresa: Eiffestrasse 80, 20537 Hamburg, Nemecko
---	---

 Distribúcia a servis pre SR	CELIMED s.r.o. Pri majeri 22, 831 06 Bratislava pozáručný servis: 02 4468 1247 e-mail: info@celimed.sk, www.celimed.sk
---	---

Please read the guide carefully before use and keep it well.
For American please refer to "°F", for European please refer to "°C".

A handheld infrared thermometer with a white body and a grey sensor housing. The device features a digital LCD screen displaying '98.6°F'. Below the screen are three small circular buttons and a larger button with a power symbol. The device has a ergonomic handle and a trigger button on the side.

User's Manual

Thank you for choosing our Infrared Thermometer. The Infrared Thermometer is used to measure an object's temperature based on the relationship between temperature and measurable infrared radiation. Simply aim the unit's probe toward the surface to be measured to obtain a quick and accurate temperature. When measuring human body, the surface temperature of the forehead and convert it to the actual human temperature by an algorithm (adjustment mode).

To ensure proper use, please be sure to read this user manual carefully. Also, please pay close attention to the safety precautions.

- In order to use this product correctly, please read the user manual before use.
- In order to properly use this product, please carefully read the full text of this manual before using, in particular the "safety precautions" section.
- Please keep the manual for easy reference.

All objects with a temperature above absolute zero emit a certain amount of infrared radiation energy based on the object's temperature. The amount of radiation energy is closely related to the wavelength distribution. When a person's forehead temperature is 96.8-98.6°F (36-37°C) Fahrenheit, it emits a wavelength of 9-13 μm of infrared radiation. Based on this principle, we are able to measure a person's actual body temperature just by measuring the forehead's surface temperature.

Warning

- Use of this thermometer is not intended as a substitute for consultation with your physician. It is dangerous for user to perform a self-evaluation and self-treatment based on the measuring result. Be sure to follow doctor's instruction.
- Keep the thermometer out of the reach of children. If the battery or any other component is swallowed, immediately call a poison center or doctor.
- Do not dispose of the battery in fire. Must be recycled or disposed properly.

- The device is precision instrument, don't drop, tramp or impose any vibration or impact on the thermometer.
- Do not touch the lens of the probe with your fingers, and do not disassemble the device by yourself.
- Before measuring temperature, please make sure the area is clear of hair or moisture (i.e. sweat).
- If you plan on measuring your temperature after exercising, eating, or bathing, please rest in a room temperature area for about 30 minutes prior to measuring the temperature for accurate results.
- For reliable and stable measurement data, please have the device placed in a room temperature area for at least 30 minutes prior to using.
- To measure temperature continuously over a duration of time, please make sure to have at least 1 minute in-

The diagram illustrates the CT-5000 digital thermometer from two perspectives: front and back. The front view on the left shows the measuring probe at the top and the Power/Scan button on the handle. The back view on the right shows the LCD display, Tone selection button, C°/F switch, Environment/Body temperature switch button, and Memory button. A large 'CELL' watermark is visible in the bottom right corner of the page.

Measuring probe

Power/Scan button

LCD

Tone selection button

C°/F switch

Environment/Body temperature switch button

Memory button

Solution	Dilution with water	Method
Ethanol	70%~80%	Wipe the probe enclosure and button with the soaked cloth for 30s.

Isopropyl Alcohol	70%~90%	As above.
Glutaraldehyde	0.5%~2%	As above.

11. Calibration

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instruction, periodic re-adjustment is not required. If any time your question the accuracy of measurement, please contact distributor or manufacturer, the contact information see last page.

12. Storage

- 
- 1) Don't put the thermometer under direct sunlight, high temperature, moist environments, flammable areas, or areas susceptible to vibration.


2) Remove the batteries if the device is not to be used for long durations.

13. Accessories

Only use original accessories. Check that the contents of the delivery are complete.

Quantity	Parts
1pc	Infrared Thermometer
2pcs	AAA batteries
1pc	User Manual

14. Trouble-Shooting

Troubles or error message	Checklists or situation	Countermeasures or solution
No response/ Automatically reset	The batteries are used up?	Replace new batteries.
	Battery in wrong polarity or type?	Take out the batteries and replace new ones.
	Poor battery contact	Take out batteries and reinsert them in the correct orientation
The thermometer show the symbol “Hi”	Temperature hampered by an air flux.	Please leave the area and wait for 30 minutes to measure. Re-measure
	In the forehead measurement mode: --Temperature readings too close together. --Measured the other object, such as the sunlight, the air from the fireplace. Hi: Higher than 109.2°F (42.9°C);	according to the manual.
	In the object measurement mode: --Temperature readings too close together. --The object temperature is higher than 212°F (100.0°C). Hi: Higher than 212°F (100.0°C);	
The thermometer show the symbol “Lo”	The hair and sweat prevent the temperature achievement.	according to the manual.
	In the forehead measurement mode: --The measuring distance is too far. --Measured the other object, such as the air from the air conditioner. Lo: Less than 89.6°F (32.0°C)	
	In the object measurement mode: --The measuring distance is too far. --Have water vapor condenses on the lens. Lo: Less than 32°F (0°C)	
	The ambient temperature is beyond of range of measurement (10°C~40°C or 50°F~104°F)	Keep the thermometer in the room Whose Temperature is (10°C~40°C or 50°F~104°F) for 30 minutes
	The sensor or hardware is damaged	Excluding the possibility of temperature allowance first, then send the device to your dealer for repair
	Low batter, but you can still use it	Keep an eye on power and continue to use.
	Lower battery, however you can't use it	Replace the new battery.

15. Specifications

Device name	Infrared Thermometer
Model	FDIR-V22
Measurement mode	Forehead and object temperature measurement modes
Measurement site	Forehead temperature
Measuring distance	1cm ~ 6cm (0.4in ~ 2.5in)
Powersupply	d.c.3V, 2×1.5V AAA batteries
Measuring range	For forehead temperature: 89.6°F ~ 109.2°F (32.0°C ~ 42.9°C) For object surface temperature: 32°F ~ 212°F(0°C ~100°C)
Measuring accuracy	For forehead temperature: 35.0°C ~ 42.0°C ± 0.2°C (95.0°F ~ 107.6°F ± 0.4°F) Outside this range: ± 0.3°C (±0.5°F)
Clinical repeatability	Within ±0.5°F (±0.3°C)
Resolution of display	0.1°F/0.1°C
Operation condition temperature	Temperature: 50.0 °F ~ 104.0 °F (10.0 °C ~ 40.0 °C) Relative humidity: ≤ 95%RH Atmospheric pressure: 70 kPa ~ 106 kPa
Transport/Storage condition temperature	Temperature: -13.0 °F ~ 131.0 °F (-25.0 °C ~ 55.0 °C) Relative humidity: ≤ 95%RH, non-condensing Atmospheric pressure: 70 kPa ~ 106 kPa
Size	149x85x44mm
Weight	110g
Expected lifetime	5 years
High body temperature hint	≥ 100.4°F (≥38°C)
Grade of waterproof	IP22
Electric shock	Internally powered ME equipment
Applied part	Type BF applied part, including the whole unit
Mode of operation	Continuous operation
Software Version	V1.0

***The above specifications are subject to change without prior notice.**
Note: ASTM laboratory accuracy requirements in the display range of 96.8°F to 102.2°F(36°C to 39°C)for this thermometer is ±0.4°F(±0.2°C), whereas for mercury-in-glass thermometer, the requirement per ASTM standards E 667-86 is ±0.2°F(±0.1°C).

16. Disposal

- 
- At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.
Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment.
Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

NOTES:

- Handing of battery and wastes method, please act according to the native law to proceed to handle.
- Take out the battery if you are not going to use the unit for a long time.

To protect the environment, dispose of empty battery at your retail store or at appropriate collection sites according to national or local regulations.

17. Normalized Symbols

 1639	Complies with the European Medical Device Regulation (EU) 2017/745, Notified Body is SGS Belgium NV.
	Authorized representative in the European Community.
	Follow the operating instructions.
	General warning sign
	Caution
	Type BF applied parts
	Batch code
IP22	IP code of the device: this device's grade of against ingress of solid foreign objects
	Disposal in accordance with Directive 2002/96/EC (WEEE)
	Medical Device
	Manufacturer information:
	Unique device identifier
	Date of manufacture

18. Electromagnetic Compatibility Information

Guidance and manufacturer’s declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic emission		
The FDIR-V22 is intended for use in the electromagnetic environment specified below. The customer or the user of FDIR-V22 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The FDIR-V22 uses RF energy only for its internal function. There for, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11 Harmonic emissions IEC 61000-3-2	Class B Class A	The FDIR-V22 suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration – electromagnetic immunity			
The FDIR-V22 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V22 should assure that it is used in such an environment.			
Immunity test	IEC 60601	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U _i ; 0,5 cycle U _i At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _i ; 1 cycle and 70 % U _i ; 25/30 cycles Single phase: at 0° 0 % U _i ; 250/300 cycle	0 % U _i ; 0,5 cycle U _i At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U _i ; 1 cycle and 70 % U _i ; 25/30 cycles Single phase: at 0° 0 % U _i ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FDIR-V22 requires continued operation during power mains interruptions, it is recommended that the FDIR-V22 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note :U _i is the a. c. mains voltage prior to application of the test level.			

Guidance and manufacturer’s declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM

Guidance and manufacturer’s declaration – electromagnetic immunity			
The FDIR-V22 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V22 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	3V 150 kHz to 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the FDIR-V22 , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=[\frac{3.5}{V_t}] \sqrt{P}$ $d=[\frac{12}{V_t}] \sqrt{P}$ $d=[\frac{3.5}{E_t}] \sqrt{P}$ 80 MHz to 800 MHz $d=[\frac{12}{E_t}] \sqrt{P}$ 800 MHz to 2.7 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^b should be less than the compliance level in each frequency range. ^c Interference may occur in the vicinity of equipment marked with the following symbol: 
	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz	
Radiated RF IEC 61000-4-3	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	

Note 1： At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2： These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from struct-ures, objects and people.

A The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz. .

B Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FDIR-V22 is used exceeds the applicable RF compliance level above, the FDIR-V22 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re–orienting or relocating the FDIR-V22 .

C Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTE

Recommended separation distances between portable and mobile RF communications equipment and the FDIR-V22				
The FDIR-V22 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FDIR-V22 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FDIR-V22 as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM and amateur radio bands $d=[\frac{3.5}{V_t}] \sqrt{P}$	150 kHz to 80 MHz in ISM and amateur radio bands $d=[\frac{12}{V_t}] \sqrt{P}$	80 MHz to 800 MHz $d=[\frac{3.5}{E_t}] \sqrt{P}$	800 MHz to 2.7 GHz $d=[\frac{7}{E_t}] \sqrt{P}$
0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	35	70
For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1： At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2： These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

19. Standard List

FDIR-V22 complies with the following standards

EN ISO 15223-1	Medical device –symbols to be used with medical device labels,labeling and information to be supplied –Part 1; General requirements
EN 1041	Information supplied by the manufacturer with medical devices
IEC 60601-1	Medical electrical equipment Part 1: General requirements for basic safety and essential performance
IEC 60601-1-11	Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home equipment and medical electrical systems used in the home
IEC 60601-1-2	Medical electrical equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
ASTM E 1965	Standard Specification for Infrared Thermometer for Intermittent Determination of Patient Temperature
IEC 62304	Medical device software - Software life-cycle processes
IEC 62366	Medical devices — Application of usability engineering to medical devices IEC 62366:2007
ISO 80601-2-56	Medical electrical equipment – Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

20.Warranty

In case of any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state.

We provide one year warranty starting from the date of purchase. Please refer to the following situations that are excluded from the free repair services within the warranty period.

- All damages caused by disassembly and repair of the device by yourselves.
 - All damages caused by dropping the device during usage, or transport.
 - All damages caused by improper usage of the device and not following the instructions on the user manual.
- Please contact after-sales service and support and enclose your product purchase receipt while claiming for warranty services.

Location of purchase:

Contact number:

Date of purchase:

Famidoc Technology Co., Ltd.
Room 401, 501, 601, 901, Building 6, Zhongtang Tian'an Cyber Park, No. 89, Jinyuan Road, Zhongtang Town, Dongguan City, Guangdong Province, 523220, P.R. China
Tel.: +86-769-89272488
Fax: +86-769-89272498
Website: www.famidoc.com

Name: Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany