

Small but Smart

Hawk-IA1 Infusion Alarm



Easy for use:

Simple operation and installation within a clip and compact light weight design allows easy movement.



Intelligent:

Compatible with various liquids and IV sets.



Prominent Precision

Accurate detection and synchronous warning help to ease patient and caregivers' burden and improve patients' safety.



Excellent battery performance:

Durable battery with long standby time saves caregivers' time and ease them from frequent recharge task.



Superior material :

Medical grade material and strong anti-interference capability assures the reliable alarm without being affected by external disturbance.



Waterproof:

Highly Waterproof and moisture-proof assures it to be used in hospital and clinic safely.

TECHNICAL SPECIFICATIONS

• **Input Power:** $\leq 2\text{VA}$

• **Alarm:** Bubble alarm, low battery alarm

• **Output Voltage:**
DC 5V

• **Alarm Volume:** $\geq 65\text{dB (A)}$

• **Input Voltage:**
AC 100V–240V 50Hz/60Hz

• **Waterproof:** IP24

• **Bubble Alarm:** $\leq 20\text{mm}$

• **Working Temperature:** 5°C-40°C

• **Weight:** 32.8g

• **Working humidity:** 10%-90%

• **Size:** 56*42*22mm

• **Battery Life:**
40 hours continuous monitoring;
4 hours continuous bubble alarm.



🏢 Shenzhen hawk medical Instrument Co., Ltd.

☎ 0755-83151901 ✉ szhk@hawkmedical.cn

🌐 www.hawkmedical.cn

📍 Add: Building C, Jianyetai Industrial Zone, No.11 Minhuan Road,
Fukang Community, Longhua Street, Shenzhen, GD, P.R.China.